



Spare sponge for CristalProfi m internal filter







- Clean, healthy water: spare sponge for the JBL CristalProfi m Internal Filter
- Easy to use: simply press the sponge into the internal filter from the outside, fully washable.
- Safe: fine pores small fish and shrimps cannot be sucked in.
- Pore size of the foam: 35 ppi
- Package contents: spare sponge for CristalProfi m Internal Filter, CristalProfi m greenline FilterPad



You may also be interested in

You can find a complete overview here: https://www.jbl.de/qr/60967





JBL CristalProfi m greenline Module Filter module to extend CristalProfi m









## **Product** information

### Safe:

The fine pores of the sponge protect small fish and shrimps from being sucked into the filter.

Just press the spare sponge into the JBL CristalProfi m Internal Filter.

### Quality of foam:

Foams used in aquaristics with live animals must always be subjected to a pollutant test, as most plastics release problematic plasticisers into the water. All the plastics used by JBL have passed this test and do NOT release any harmful substances into the water.

## Cleaning of plastics:

For many months it is sufficient to rinse the dirty filter sponges under lukewarm running water and then inoculate them with filter bacteria (JBL FilterStart) before using them again. However, as mechanical cleaning NEVER leads to complete cleaning and dirt will remain in the filter sponge, it is recommended you replace the filter sponges about every 6 months.

If you place a freshly cleaned filter sponge in a glass of clean water overnight, you will be able to smell how much dirt is still present in the supposedly clean sponge the next morning!

Inoculating filter material with cleansing bacteria: Initial colonisation of filter material with cleansing bacteria takes about 1-2 weeks in freshwater and twice as long in seawater. If the filter material was already in use, the time for a "recolonisation" is doubled, as recolonisation does not take place so easily on dead bacterial lawns. The addition of concentrated cleaning bacteria such as in JBL FilterStart shortens the biological start-up time to about 24 hours.

Another reason is to fill the existing settlement space in the filter material with beneficial bacteria so that pathogenic germs find no place to spread. This has now been scientifically proven.

Further information	
FAQ	~
Blog	~
Press	~
Laboratory/calculator	×
Worth reading	~
Spare parts	×
Video	~
GarantiePlus	×
Instructions	~
QR code	



30.12.2023 Date: Produced by:





# JBL CP m greenline FilterPad Product details



Article data	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
EAN number	4014162609670
EAN as barcode	
Content	2 piece(s)
-	-
Expiry months	-
RRP incl. VAT	8,97 €
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	38 g
Net weight	10 g
Weight change	1000

Disposal	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	-
PPK	42 g
Plastic small	-
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700



30.12.2023 Date: Produced by:

# TECHNICAL PRODUCTS TECHNIK TECHNIQUE





Features	
Animal species	Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichilds (South America), Clawed frogs, Clawed shrimps, Corals, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle
Animal size	1 - 5 cm
Animal age group	Small aquarium fish
Volume of habitat	20 - 80 (20 - 100 cm)
Material	polyether
Food type	-
Colour	black
Dosage	-
Transport conditions	-







Electronic label / illuminant	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-
Hose connections suction/in	-
Delivery head	-



Date: 30.12.2023 Produced by:





Food type	-
Sub product type	-
Dosing	-



Date: 30.12.2023 Produced by:







## Additional information for the specialist trade sector

Article data	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
VAT	19%
Sales unit (SU)	3
Volume packaging	0.481
Dimensions (I/w/h)	35 mm/100 mm/137 mm
Layer	270
Pallet	540
Category of products	3
Customs tariff	84219990
Country of origin	CN
Type of packaging	Faltschachtel/Karton

PU 1 data	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
PU 1 material	film gr
PU 1 weight	2 g
PU 1 lengh	150 mm
PU 1 width/depth	135 mm
PU 1 height	140 mm

PU 2 data	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
PU 2 material	cardboard gr
PU 2 weight	1050 g
PU 2 lengh	730 mm
PU 2 width/depth	320 mm
PU 2 height	440 mm

Trade data	
Product name	JBL CP m greenline FilterPad
Art. No.	6096700
Till receipt text	Profi m Pad
Shelf placement	-



Date: 30.12.2023 Produced by: