Quick test to determine the acidity

Suitable for: < 🖈 🏏 🛸

CARE

PFLEGE ENTRETIEN



- Easy and reliable monitoring of the water values of aquariums and ponds. Determines the ideal pH values for freshwater and seawater
- Quick test: fill plastic vial with sample water, add reagent, read value of the sample off colour chart
- When to use: for setting up a new aquarium: daily for one week, afterwards weekly
- Online Laboratory: regular control for a healthy aquarium/pond with natural conditions. JBL has water tests in the form of quick tests or colour change tests for every water analysis. Testing aquarium water for healthy, clear water
- Contents: 1 refill pack for quick test, contents suffice for approx. 50 measurements, 1 reagent pH 3.0-10.0



You may also be interested in You can find a co

You can find a complete overview here: https://www.jbl.de/qr/11111



JBL Test Combi Set plus Fe Test case with most important water tests + iron test



JBL Testlab Test case with 13 tests f. freshwater analysis Testlab



JBL Testlab Marin Professional test case for the analysis of saltwater



JBL EasyTest 6in1 Test strips for quick aquarium water testing











JBL component set for water tests









JBL Testlab empty + inserts





Product information

Healthy aquarium/healthy pond with conditions close to nature The right water values are dependent on the fish stock and the plants in the aquarium/pond. Even if the water looks clear it can be contaminated. With bad values diseases or algae can appear in the aquarium/pond. To maintain a healthy aquarium/pond with conditions close to nature it is important to check and adapt the water values regularly.

For each water analysis JBL provides water tests in the form of quick tests or colour change tests. These determine a certain value or several values in one go. With these water tests you can recognise algae problems and negative nitrate, nitrite, potassium, magnesium values etc.

Why test?

The well-being of aquarium and pond dwellers and the growth of aquatic plants depend to a large extent on the pH level being kept as constant as possible. The pH level can also affect many substances dissolved in the water. Fluctuations in the pH level are definitely to be avoided.

JBL Online Laboratory

Regular control for a healthy aquarium/pond with conditions close to nature. Enter your water values into the JBL Online Laboratory and get an in-depth analysis of your values within seconds.

Recommended pH values: Freshwater aquarium (community aquarium): 6.5-7.5 Lake Malawi/Lake Tanganyika aquarium: 7.8-9.2 Plant aquarium with few fish (aquascaping): 6.0-7.0 Marine aquarium: 7.9-8.5 Pond: 7.5-8.5

| Further information | |
|-----------------------|---|
| FAQ | ✓ |
| Blog | ~ |
| Press | ~ |
| Laboratory/calculator | × |
| Worth reading | ✓ |
| Spare parts | ~ |
| Video | ~ |
| GarantiePlus | × |
| Instructions | ~ |
| QR code | |



Article data Product name JBL pH Test Set 3.0-10.0 JBL pH Reagent 3.0-10.0 Art. No. 2534200 2534300 EAN number 4014162253422 4014162253439 EAN as barcode For 50 Tests Expiry months 48 48 RRP incl. VAT Base price Nominal filling quantity Base quantity 1 1 Gross weight 42 g 27 g Net weight 27 g 10 g 1000 1000 Weight change Die . 1

| Disposal | | |
|-------------------------|--------------------------|-------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| Green dot | ✓ | ✓ |
| Group electronic waste | - | - |
| Disposal weight | - | - |
| Battery type | - | - |
| Battery return | - | - |
| Battery rechargeable | - | - |
| Disposal weight battery | - | - |
| Non-returnable glass | - | - |
| РРК | 14 g | 7 g |
| Plastic small | 6 g | - |
| Plastic large | - | - |
| Disposal weight metal | 0 g | 0 g |

| Features | | |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| Animal species | Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle | Arowana, Axoloti, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle |





| Features | | |
|----------------------|-----------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| Animal size | For all animal sizes | For all animal sizes |
| Animal age group | All aquarium fish | All aquarium fish |
| Volume of habitat | 80 tests | 80 tests |
| Material | Methyl red, bromo xylenol blue, o- cresolphthalein complexion, phenolphthalein ethanol, purified water | Methyl red, bromo xylenol blue, o- cresolphthalein complexion, phenolphthalein ethanol, purified water |
| Food type | - | - |
| Colour | red | red |
| Dosage | 5 ml sample water, add 4 drops reagent, mix briefly and let stand for 3 minutes | 5 ml sample water, add 4 drops reagent, mix briefly and let stand for 3 minutes |
| Transport conditions | - | - |



Date: 30.12.2023 Produced by:



| Electronic label / illur | ninant | |
|--------------------------|--------------------------|-------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| Ambient temperature | - | - |
| Start time | - | - |
| Mercury | - | - |
| Tube length | - | - |
| Service life | - | - |
| Lumen | - | - |
| CRI value | - | - |
| Dimmable | - | - |
| Switching cycles | - | - |
| PAR value | - | - |
| Energy efficiency class | - | - |
| UV-A | - | - |
| UV-B | - | - |
| UV-C | - | - |
| Colour temperature | - | - |
| Base designation | - | - |
| | | |
| Technical data | | |
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| Range in litres | - | - |
| Range from - to | - | - |
| Range in days | - | - |
| Range tank length | - | - |
| Output in watts | - | - |
| Output per hour | - | - |
| Output per day | - | - |
| Height | - | - |
| Length | - | - |
| Width | - | - |
| - | | |

| Output per day | - | - |
|-------------------------------|---|---|
| Height | - | - |
| Length | - | - |
| Width | - | - |
| Diameter | - | - |
| Voltage | - | - |
| For | - | - |
| T8 26mm (watt) | - | - |
| T5 16mm (watt) | - | - |
| Size | - | - |
| Content for | - | - |
| Filter container volume | - | - |
| Volume filter media | - | - |
| Hose connections pressure/out | - | - |
| Hose connections suction/in | - | - |
| Delivery head | - | - |
| | | |







| Food type | |
|------------------|---------------------------------------------------------------------------------|
| Sub product type | |
| Dosing | 5 ml sample water, add 4 drops reagent, mix briefly and let stand for 3 minutes |





| - safety information in accordance with GHS | | |
|---------------------------------------------|------------------------------------------------------------------------------------------------------|--|
| Signal word | Warning | |
| Hazard symbol | | |
| Risk phrase | H226: Flammable liquid and vapour. | |
| safety note | P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. | |

| Further safety information | | |
|----------------------------|--------------|--|
| Safety instructions | \checkmark | |
| QR code | | |





Additional information for the specialist trade sector

| Article data | | |
|----------------------|--------------------------|-------------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| VAT | 19% | 19% |
| Sales unit (SU) | 6 | 1 |
| Volume packaging | 0.31 | 0.451 |
| Dimensions (l/w/h) | 35 mm/70 mm/125 mm | 30 mm/107 mm/140 mm |
| Layer | 960 | 0 |
| Pallet | 2880 | 0 |
| Category of products | 1 | 1 |
| Customs tariff | 38220000 | 38220000 |
| Country of origin | DE | DE |
| Type of packaging | Faltschachtel/Karton | Skin / blister card packaging |

| PU 1 data | | |
|------------------|--------------------------|-------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| PU 1 material | film gr | cardboard gr |
| PU 1 weight | 4.3 g | 660 g |
| PU 1 lengh | 215 mm | 580 mm |
| PU 1 width/depth | 35 mm | 390 mm |
| PU 1 height | 185 mm | 310 mm |

| PU 2 data | | |
|------------------|--------------------------|-------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| PU 2 material | cardboard gr | None |
| PU 2 weight | 500 g | - |
| PU 2 lengh | 572 mm | - |
| PU 2 width/depth | 377 mm | - |
| PU 2 height | 410 mm | - |

| Trade data | | |
|-------------------|--------------------------|-------------------------|
| Product name | JBL pH Test Set 3.0-10.0 | JBL pH Reagent 3.0-10.0 |
| Art. No. | 2534200 | 2534300 |
| Till receipt text | pH Test 3-10 | pH Test 3-10 Re |
| Shelf placement | - | - |

