



# JBL Testlab Marin

Professional test case for the analysis of saltwater

Suitable for:



- Simple and reliable monitoring of the water values in marine aquariums. Determination of: pH value, carbonate hardness, calcium, magnesium, copper, ammonium, nitrite, nitrate, phosphate, silicate, oxygen
- Easy to use: syringe for exact measurement of the water quantity, detailed instructions for use
- Laboratory comparator system to compensate the inherent water colouring, analysis record sheet, instructions for problem solving
- Child-proof reagent bottles, water-proof plastic case, refillable kit
- Package contents: 1 professional test case for analysing seawater. Incl. 11 reagents, 6 glass vials, 2 syringes, thermometer, ball pen, analysis record sheet, laboratory comparator system, detailed instructions for use

You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/25503>



**JBL pH Test 3.0-10.0**  
Quick test to determine the acidity



**JBL pH Test 6.0-7.6**  
pH test for freshwater aquariums in the range 6.0-7.6



**JBL pH Test 7.4-9.0**  
pH test in aquariums and ponds in the range 7.4-9.0



**JBL GH Test**  
Quick test to determine the total hardness





# JBL Testlab Marin

## Accessories



**JBL Aquarium  
Thermometer Float**  
Aquarium thermometer



# JBL Testlab Marin

## Spare parts



**JBL assortment Test Colour Charts**  
Colour charts for JBL water tests



**JBL component set for water tests**



**JBL Testlab empty + inserts**



# JBL Testlab Marin

## Product information

Healthy aquariums with conditions close to nature

The right water values are dependent on the fish stock, and the algae and corals in the aquarium. Even if the water looks clear it can be contaminated. With bad values diseases or algae can appear in the aquarium. To maintain a healthy aquarium with conditions close to nature it is important to check and adapt the water values regularly. With the JBL Testlab Marin you can determine the most important water values in your aquarium.

Precise measurement of the following values:

pH test: acidity of the water from 7.4 to 9.0

KH test: pH stability of the water (carbonate hardness)

Ca test: most important macro-element for the growth of invertebrates

Mg test: important macro-element for the growth of invertebrates and calciferous red algae

Cu test: deadly heavy metal for invertebrates. Also important for setting the therapeutic dose of medications

NO<sub>2</sub> test: toxic nitrogen compound

NO<sub>3</sub> test: too high nitrate values cause unwanted algae growth


NH<sub>4</sub> test: indication of non-toxic ammonium, determination of the toxic ammonia using table

PO<sub>4</sub> test: high phosphate values lead to unwanted algae growth

SiO<sub>2</sub> test: cause for diatoms

O<sub>2</sub> test: basic requirement for fish and plants

You can analyse the water values you have measured in the JBL Online Laboratory free of charge and receive recommendations.

Further information	
FAQ	✓
Blog	✓
Press	✓
Laboratory/calculator	✗
Worth reading	✓
Spare parts	✓
Video	✗
GarantiePlus	✗
Instructions	✓
QR code	





# JBL Testlab Marin



## Product details

Article data	
Product name	JBL Testlab Marin
Art. No.	2550300
EAN number	4014162255037
EAN as barcode	
Content	11 test kits
-	-
Expiry months	24
RRP incl. VAT	-
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	1300 g
Net weight	1300 g
Weight change	1000

Disposal	
Product name	JBL Testlab Marin
Art. No.	2550300
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	39 g
PPK	-
Plastic small	256 g
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL Testlab Marin
Art. No.	2550300
Animal species	Corals, Gobies, Mussels, Snails
Animal size	For all animal sizes
Animal age group	All aquarium fish
Volume of habitat	For all aquariums
Material	-
Food type	-





Features	
Colour	multicoloured
Dosage	see instruction manual
Transport conditions	-





# JBL Testlab Marin

Electronic label / illuminant	
Product name	JBL Testlab Marin
Art. No.	2550300
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL Testlab Marin
Art. No.	2550300
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-
Hose connections suction/in	-
Delivery head	-





# JBL Testlab Marin

Food type	-
Sub product type	-
Dosing	see instruction manual







# JBL Testlab Marin

## Safety information

- safety information in accordance with GHS	
Signal word	Danger
Hazard symbol	
Risk phrase	<p>H261: In contact with water releases flammable gases.                      EUH206: Warning! Do not use together with other products. May release dangerous gases (chlorine).                      EUH032: Contact with acids liberates very toxic gas.                      H412: Harmful to aquatic life with long lasting effects.                      H400: Very toxic to aquatic life.                      H335: May cause respiratory irritation.                      H341: Suspected of causing genetic defects.                      H317: May cause an allergic skin reaction.                      H314: Causes severe skin burns and eye damage.                      H225: Highly flammable liquid and vapour.</p>
safety note	<p>P101: If medical advice is needed, have product container or label at hand.                      P102: Keep out of reach of children.                      P103: Read label before use.                      P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.                      P301+P310a: IF SWALLOWED: Immediately call a POISON CENTER/ doctor.                      P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.                      P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.                      P405: Store locked up.                      P501a: Dispose of contents/container in accordance with local/regional/national/international regulations.</p>

Further safety information	
Safety instructions	✓
QR code	





## Additional information for the specialist trade sector

Article data	
Product name	JBL Testlab Marin
Art. No.	2550300
VAT	19%
Sales unit (SU)	1
Volume packaging	10.32l
Dimensions (l/w/h)	85 mm/382 mm/318 mm
Layer	24
Pallet	72
Category of products	1
Customs tariff	38220000
Country of origin	DE
Type of packaging	case

PU 1 data	
Product name	JBL Testlab Marin
Art. No.	2550300
PU 1 material	cardboard gr
PU 1 weight	1700 g
PU 1 length	700 mm
PU 1 width/depth	410 mm
PU 1 height	370 mm

PU 2 data	
Product name	JBL Testlab Marin
Art. No.	2550300
PU 2 material	None
PU 2 weight	-
PU 2 length	-
PU 2 width/depth	-
PU 2 height	-

Trade data	
Product name	JBL Testlab Marin
Art. No.	2550300
Till receipt text	Testlab Marin
Shelf placement	-