



# JBL PROCRISTAL i30 FilterSponge

Replacement filter sponge for CristalProfi i 30

Suitable for: 



- Clear water, healthy fish: replacement filter sponge for JBL ProCristal i 30
- Clear aquarium water: replace the filter sponge every 3 months.
- Best cleaning: pore size 30 ppi
- Best biological filtering: promotes the settlement of beneficial cleansing bacteria.
- Package contents: replacement sponge for JBL CristalProfi i 30, ProCristal i30 FilterSponge



You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/60993>



## JBL PROCRISTAL i30 FilterModul

Filter module for ProCristal i30





# JBL PROCRISTAL i30 FilterSponge



## Product information

### Clear aquarium water

Replacing the filter sponge regularly ensures clear aquarium water. The water is perfectly cleaned thanks to the pore size of 30 ppi. At the same time the sponge promotes the settlement of beneficial cleansing bacteria for an ideal biological filtering.

### Quality of foam:

Foams used in aquaristics with live animals must always be subjected to a pollutant test, as most plastics release problematic plasticisers into the water. All the plastics used by JBL have passed this test and do NOT release any harmful substances into the water.


### Cleaning of plastics:

For many months it is sufficient to rinse the dirty filter sponges under lukewarm running water and then inoculate them with filter bacteria (JBL FilterStart) before using them again. However, as mechanical cleaning NEVER leads to complete cleaning and dirt will remain in the filter sponge, it is recommended you replace the filter sponges about every 6 months.

### Inoculating filter material with cleansing bacteria:

Initial colonisation of filter material with cleansing bacteria takes about 1-2 weeks in freshwater and twice as long in seawater. If the filter material was already in use, the time for a "recolonisation" is doubled, as recolonisation does not take place so easily on dead bacterial lawns. The addition of concentrated cleaning bacteria such as in JBL FilterStart shortens the biological start-up time to about 24 hours.

Another reason is to fill the existing settlement space in the filter material with beneficial bacteria so that pathogenic germs find no place to spread. This has now been scientifically proven.

Further information	
FAQ	✓
Blog	✓
Press	✓
Laboratory/calculator	✗
Worth reading	✓
Spare parts	✗
Video	✓
GarantiePlus	✗
Instructions	✓
QR code	





# JBL PROCRISTAL i30 FilterSponge



## Product details

Article data	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
EAN number	4014162609939
EAN as barcode	
For	PC i30
-	-
Expiry months	-
RRP incl. VAT	5,03 €
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	16 g
Net weight	1 g
Weight change	1000

Disposal	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	-
PPK	6 g
Plastic small	8.5 g
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
Animal species	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle





Features	
Animal size	1 - 5 cm
Animal age group	Small aquarium fish
Volume of habitat	10 - 40 L (15 - 40 cm)
Material	polyether
Food type	-
Colour	blue
Dosage	-
Transport conditions	-



# JBL PROCRISTAL i30 FilterSponge

Electronic label / illuminant	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-
Hose connections suction/in	-
Delivery head	-





# JBL PROCRISTAL i30 FilterSponge

Food type	-
Sub product type	-
Dosing	-





## Additional information for the specialist trade sector

Article data	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
VAT	20%
Sales unit (SU)	4
Volume packaging	0.63l
Dimensions (l/w/h)	50 mm/96 mm/140 mm
Layer	0
Pallet	0
Category of products	3
Customs tariff	84219990
Country of origin	IT
Type of packaging	Skin / blister card packaging

PU 1 data	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
PU 1 material	cardboard gr
PU 1 weight	44 g
PU 1 length	120 mm
PU 1 width/depth	100 mm
PU 1 height	145 mm

PU 2 data	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
PU 2 material	cardboard gr
PU 2 weight	642 g
PU 2 length	521 mm
PU 2 width/depth	379 mm
PU 2 height	316 mm

Trade data	
Product name	JBL PROCRISTAL i30 FilterSponge
Art. No.	6099300
Till receipt text	CP i30Sponge
Shelf placement	-

