



JBL UniBloc CristalProfi i60/80/100/200

Foam cartridge for Cristal Profi i

Suitable for:



- Clear and healthy water: replacement foam cartridge for aquarium filters CristalProfi i60-200
- Use: add FilterStart, insert into filter module
- Foam with ideal service life thanks to 20 ppi pores (ppi = pores per inch/2.5 cm), effective mechanical and biological filtering.
- No release of harmful substances: 100% water neutral synthetic material
- Package contents: 1 replacement foam cartridge (without internal strainer tube)

You may also be interested in

You can find a complete overview here: <https://www.jbl.de/qr/60928>



JBL CRISTALPROFI i greenline module

Filter module for internal filter series CristalProfi i





JBL UniBloc CristalProfi i60/80/100/200



Product information

Quality of foam:

Foams used in aquaristics with live animals must always be subjected to a pollutant test, as most plastics release problematic plasticisers into the water. All the plastics used by JBL have passed this test and do NOT release any harmful substances into the water.

Cleaning of plastics:


For many months it is sufficient to rinse the dirty filter sponges under lukewarm running water and then inoculate them with filter bacteria (JBL FilterStart) before using them again. However, as mechanical cleaning NEVER leads to complete cleaning and dirt will remain in the filter sponge, it is recommended you replace the filter sponges about every 6 months.

If you place a freshly cleaned filter sponge in a glass of clean water overnight, you will be able to smell how much dirt is still present in the supposedly clean sponge the next morning!

Inoculating filter material with cleansing bacteria:

Initial colonisation of filter material with cleansing bacteria takes about 1-2 weeks in freshwater and twice as long in seawater. If the filter material was already in use, the time for a "recolonisation" is doubled, as recolonisation does not take place so easily on dead bacterial lawns. The addition of concentrated cleaning bacteria such as in JBL FilterStart shortens the biological start-up time to about 24 hours.

Another reason is to fill the existing settlement space in the filter material with beneficial bacteria so that pathogenic germs find no place to spread. This has now been scientifically proven.

Further information	
FAQ	✓
Blog	✓
Press	✗
Laboratory/calculator	✗
Worth reading	✓
Spare parts	✗
Video	✓
GarantiePlus	✗
Instructions	✓
QR code	





JBL UniBloc CristalProfi i60/80/100/200

Product details

Article data	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
EAN number	4014162609281
EAN as barcode	
Content	1 piece(s)
For CPi	60-200
Expiry months	-
RRP incl. VAT	7,03 €
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	36 g
Net weight	4 g
Weight change	1000

Disposal	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	-
PPK	32 g
Plastic small	-
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800





Features	
Animal species	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Corals, Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, freshwater butterflyfish, turtle
Animal size	2 - 10 cm
Animal age group	Small and medium-sized aquarium fish
Volume of habitat	200 L
Material	polyether
Food type	-
Colour	blue
Dosage	-
Transport conditions	-



JBL UniBloc CristalProfi i60/80/100/200

Electronic label / illuminant	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-





Technical data	
Hose connections suction/in	-
Delivery head	-





JBL UniBloc CristalProfi i60/80/100/200

Food type	-
Sub product type	-
Dosing	-





Additional information for the specialist trade sector

Article data	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
VAT	24%
Sales unit (SU)	3
Volume packaging	0.55l
Dimensions (l/w/h)	74 mm/74 mm/100 mm
Layer	630
Pallet	1260
Category of products	2
Customs tariff	39269097
Country of origin	CN
Type of packaging	Faltschachtel/Karton

PU 1 data	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
PU 1 material	film gr
PU 1 weight	2.6 g
PU 1 length	75 mm
PU 1 width/depth	220 mm
PU 1 height	100 mm

PU 2 data	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
PU 2 material	cardboard gr
PU 2 weight	650 g
PU 2 length	400 mm
PU 2 width/depth	420 mm
PU 2 height	340 mm

Trade data	
Product name	JBL UniBloc CristalProfi i60/80/100/200
Art. No.	6092800
Till receipt text	UB CPI80-200
Shelf placement	-

