



JBL Acclimol

Water conditioner for acclimatisation

Suitable for: 



- For aquarium water ideally suited to the requirements of the fish: when adding new fish, after working in the aquarium, during fish transport to activate the resistance of freshwater aquarium fish
- Protection against stress and support of the immune system thanks to multivitamin complex with vitamin C
- Healthy water, healthy fish: protection of mucous membranes, gills, skin and fins by means of protective colloids
- Prevention of iodine deficiency symptoms thanks to iodine complex
- Contents: 1 bottle Acclimol. To use: 10 ml/40 l water. Example: in a 60 l aquarium: 15 ml for setting up a new aquarium, 5 ml for 1/3 water change every 2-3 weeks



JBL Acclimol

Accessories



JBL dosing pump Acclimol/AccliPond 500ml

For the fish protection during
transport





JBL Acclimol

Product information

Nature at Home – Home for Fish and Plants


Fish need to be as comfortable in the aquarium as they are in nature. For the water change, aquarium setup and to strengthen the fish after disease. Relocation causes stress. The transport and insertion into an aquarium is highly stressful for fish. With JBL's Acclimatisation Protection you can reduce the stress for the fish and ease their adaptation to the environment.

Advantages of the JBL Acclimatisation Protection

Reduces stress, increases their resistance against disease, reduces fish loss, protects the mucous membranes, gills, skin and fins thanks to aloe vera, strengthens the immune system and activates their resistance with vitamin C complex, prevents iodine deficiency symptoms by means of its iodine complex.

Recommendation for acclimatisation

Switch off aquarium lighting, attach the opened transport bag with a clip, add aquarium water in portions into the bag every 1/4 hour, put a dosage of JBL Acclimol into the aquarium before inserting the fish, use a net to move the fish from the bag into the aquarium, let the fish swim out of the net on their own. Pour the transport water down the drain and not into the aquarium. Switch on the aquarium lighting one hour later.

| Further information | |
|-----------------------|---|
| FAQ | ✓ |
| Blog | ✓ |
| Press | ✓ |
| Laboratory/calculator | ✗ |
| Worth reading | ✓ |
| Spare parts | ✗ |
| Video | ✓ |
| GarantiePlus | ✗ |
| Instructions | ✓ |
| QR code |  |



JBL Acclimol

Product details

| Article data | | | | |
|--------------------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| EAN number | 4014162230706 | 4014162230713 | 4014162230720 | 4014162230737 |
| EAN as barcode | | | | |
| Content | 50 ml | 100 ml | 250 ml | 500 ml |
| Range | 60-200 l | 400 l | 1000 l | 2000 l |
| Expiry months | - | - | - | - |
| RRP incl. VAT | 5,36 € | 8,35 € | 17,51 € | 29,18 € |
| Base price | 107.2 € | 83.5 € | 70.04 € | 58.36 € |
| Nominal filling quantity | 50 ml | 100 ml | 250 ml | 500 ml |
| Base quantity | 1 l | 1 l | 1 l | 1 l |
| Gross weight | 63 g | 124 g | 300.2 g | 578.7 g |
| Net weight | 50 g | 100 g | 252 g | 502 g |
| Weight change | - | 1000 | 1000 | 1000 |

| Disposal | | | | |
|-------------------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| Green dot | ✓ | ✓ | ✓ | ✓ |
| Group electronic waste | - | - | - | - |
| Disposal weight | - | - | - | - |
| Battery type | - | - | - | - |
| Battery return | - | - | - | - |
| Battery rechargeable | - | - | - | - |
| Disposal weight battery | - | - | - | - |
| Non-returnable glass | - | - | - | - |
| PPK | - | - | - | - |
| Plastic small | 13 g | 24 g | 48.2 g | 76.7 g |
| Plastic large | - | - | - | - |
| Disposal weight metal | 0 g | 0 g | 0 g | 0 g |

| Features | | | | |
|--------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |






| Features | | | | |
|----------------------|---|---|---|---|
| Animal species | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Tetra, Veiltails, blood parrot cichlids, freshwater butterflyfish | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Tetra, Veiltails, blood parrot cichlids, freshwater butterflyfish | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Tetra, Veiltails, blood parrot cichlids, freshwater butterflyfish | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Tetra, Veiltails, blood parrot cichlids, freshwater butterflyfish |
| Animal size | For all animal sizes | For all animal sizes | For all animal sizes | For all animal sizes |
| Animal age group | All aquarium fish | All aquarium fish | All aquarium fish | All aquarium fish |
| Volume of habitat | - | 400 L | 1,000 L | 2,000 L |
| Material | Aloe lingua exe, immunostimulans (based on alginic acid), iodine, vitamin C, purified water | Aloe lingua exe, immunostimulans (based on alginic acid), iodine, vitamin C, purified water | Aloe lingua exe, immunostimulans (based on alginic acid), iodine, vitamin C, purified water | Aloe lingua exe, immunostimulans (based on alginic acid), iodine, vitamin C, purified water |
| Food type | - | - | - | - |
| Colour | yellow | yellow | yellow | yellow |
| Dosage | Always calculate for the complete aquarium content: 10 ml per 40 l water. During transport: 1 ml per 3 l water | Always calculate for the complete aquarium content: 10 ml per 40 l water. During transport: 1 ml per 3 l water | Always calculate for the complete aquarium content: 10 ml per 40 l water. During transport: 1 ml per 3 l water | Always calculate for the complete aquarium content: 10 ml per 40 l water. During transport: 1 ml per 3 l water |
| Transport conditions | - | - | - | - |





| Article data | |
|--------------------------|---|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| EAN number | 4014162200563 |
| EAN as barcode |  |
| Content | 5000 ml |
| Range | 20000 l |
| Expiry months | - |
| RRP incl. VAT | 102,51 € |
| Base price | 20.5 € |
| Nominal filling quantity | 5 l |
| Base quantity | 1 l |
| Gross weight | 5250 g |
| Net weight | 5010 g |
| Weight change | 1000 |

| Disposal | |
|-------------------------|------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| Green dot | ✓ |
| Group electronic waste | - |
| Disposal weight | - |
| Battery type | - |
| Battery return | - |
| Battery rechargeable | - |
| Disposal weight battery | - |
| Non-returnable glass | - |
| PPK | - |
| Plastic small | 240 g |
| Plastic large | - |
| Disposal weight metal | 0 g |

| Features | |
|-------------------|---|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| Animal species | Armored catfish, Arowana, Axolotl, Barbels, Bettas, Bichirs/ reedfish, Blowfish, Catfish, Cichlids (South America), Clawed frogs, Clawed shrimps, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Tetra, Veiltails, blood parrot cichlids, freshwater butterflyfish |
| Animal size | For all animal sizes |
| Animal age group | All aquarium fish |
| Volume of habitat | 20,000 L |





| Features | |
|----------------------|--|
| Material | Aloe lingua ex, immunostimulans (based on alginic acid), iodine, vitamin C, purified water |
| Food type | - |
| Colour | yellow |
| Dosage | Always calculate for the complete aquarium content: 10 ml per 40 l water. During transport: 1 ml per 3 l water |
| Transport conditions | - |





JBL Acclimol

| Electronic label / illuminant | | | | |
|-------------------------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| Ambient temperature | - | - | - | - |
| Start time | - | - | - | - |
| Mercury | - | - | - | - |
| Tube length | - | - | - | - |
| Service life | - | - | - | - |
| Lumen | - | - | - | - |
| CRI value | - | - | - | - |
| Dimmable | - | - | - | - |
| Switching cycles | - | - | - | - |
| PAR value | - | - | - | - |
| Energy efficiency class | - | - | - | - |
| UV-A | - | - | - | - |
| UV-B | - | - | - | - |
| UV-C | - | - | - | - |
| Colour temperature | - | - | - | - |
| Base designation | - | - | - | - |

| Technical data | | | | |
|-------------------------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| Range in litres | - | - | - | - |
| Range from - to | - | - | - | - |
| Range in days | - | - | - | - |
| Range tank length | - | - | - | - |
| Output in watts | - | - | - | - |
| Output per hour | - | - | - | - |
| Output per day | - | - | - | - |
| Height | 145 mm | 125 mm | 180 mm | 220 mm |
| Length | 25 mm | 60 mm | 78 mm | 97 mm |
| Width | 25 mm | 30 mm | 33 mm | 40 mm |
| Diameter | 25 mm | - | - | - |
| Voltage | - | - | - | - |
| For | - | - | - | - |
| T8 26mm (watt) | - | - | - | - |
| T5 16mm (watt) | - | - | - | - |
| Size | - | - | - | - |
| Content for | - | - | - | - |
| Filter container volume | - | - | - | - |
| Volume filter media | - | - | - | - |
| Hose connections pressure/out | - | - | - | - |
| Hose connections suction/in | - | - | - | - |
| Delivery head | - | - | - | - |





| Electronic label / illuminant | |
|-------------------------------|------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| Ambient temperature | - |
| Start time | - |
| Mercury | - |
| Tube length | - |
| Service life | - |
| Lumen | - |
| CRI value | - |
| Dimmable | - |
| Switching cycles | - |
| PAR value | - |
| Energy efficiency class | - |
| UV-A | - |
| UV-B | - |
| UV-C | - |
| Colour temperature | - |
| Base designation | - |

| Technical data | |
|-------------------------------|------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| Range in litres | - |
| Range from - to | - |
| Range in days | - |
| Range tank length | - |
| Output in watts | - |
| Output per hour | - |
| Output per day | - |
| Height | 230 mm |
| Length | 145 mm |
| Width | 190 mm |
| Diameter | - |
| Voltage | - |
| For | - |
| T8 26mm (watt) | - |
| T5 16mm (watt) | - |
| Size | - |
| Content for | - |
| Filter container volume | - |
| Volume filter media | - |
| Hose connections pressure/out | - |
| Hose connections suction/in | - |
| Delivery head | - |





JBL Acclimol


| | |
|------------------|--|
| Food type | - |
| Sub product type | - |
| Dosing | Always calculate for the complete aquarium content: 10 ml per 40 l water. During transport: 1 ml per 3 l water |





JBL Acclimol

Safety information

| Further safety information | |
|----------------------------|--|
| Safety instructions | ✓ |
| QR code |  |





Additional information for the specialist trade sector

| Article data | | | | |
|----------------------|--------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| VAT | 25% | 25% | 25% | 25% |
| Sales unit (SU) | 24 | 6 | 6 | 3 |
| Volume packaging | 0.09l | 0.23l | 0.46l | 0.85l |
| Dimensions (l/w/h) | 25 mm/25 mm/145 mm | 30 mm/60 mm/125 mm | 33 mm/78 mm/180 mm | 40 mm/97 mm/220 mm |
| Layer | 960 | 864 | 252 | 192 |
| Pallet | 4800 | 1728 | 1260 | 768 |
| Category of products | 1 | 1 | 1 | 1 |
| Customs tariff | 38249996 | 38249996 | 38249996 | 38249996 |
| Country of origin | DE | DE | DE | DE |
| Type of packaging | Flasche | Flasche | Flasche | Flasche |

| PU 1 data | | | | |
|------------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| PU 1 material | cardboard gr | film gr | film gr | film gr |
| PU 1 weight | 85 g | 2.8 g | 4.6 g | 6 g |
| PU 1 length | 160 mm | 60 mm | 82 mm | 97 mm |
| PU 1 width/depth | 105 mm | 183 mm | 210 mm | 127 mm |
| PU 1 height | 147 mm | 125 mm | 180 mm | 220 mm |

| PU 2 data | | | | |
|------------------|-------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| PU 2 material | cardboard gr | cardboard gr | cardboard gr | cardboard gr |
| PU 2 weight | 253 g | 296 g | 332 g | 326 g |
| PU 2 length | 553 mm | 366 mm | 412 mm | 390 mm |
| PU 2 width/depth | 177 mm | 236 mm | 290 mm | 270 mm |
| PU 2 height | 165 mm | 256 mm | 186 mm | 225 mm |

| Trade data | | | | |
|-------------------|------------------------------------|---------------------|---------------------|---------------------|
| Product name | JBL Acclimol 50ml | JBL Acclimol 100 ml | JBL Acclimol 250 ml | JBL Acclimol 500 ml |
| Art. No. | 2307000 | 2307100 | 2307200 | 2307300 |
| Till receipt text | Acclimol 50ml | Acclimol 100ml | Acclimol 250ml | Acclimol 500ml |
| Shelf placement | Water conditioner, aquarium system | - | - | - |



| Article data | |
|----------------------|----------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| VAT | 25% |
| Sales unit (SU) | 1 |
| Volume packaging | 6.34l |
| Dimensions (l/w/h) | 190 mm/145 mm/230 mm |
| Layer | 32 |
| Pallet | 64 |
| Category of products | 1 |
| Customs tariff | 38249996 |
| Country of origin | DE |
| Type of packaging | Kanister |

| PU 1 data | |
|------------------|------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| PU 1 material | None |
| PU 1 weight | - |
| PU 1 length | - |
| PU 1 width/depth | - |
| PU 1 height | - |

| PU 2 data | |
|------------------|------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| PU 2 material | None |
| PU 2 weight | - |
| PU 2 length | - |
| PU 2 width/depth | - |
| PU 2 height | - |

| Trade data | |
|-------------------|------------------|
| Product name | JBL Acclimol 5 l |
| Art. No. | 2005600 |
| Till receipt text | Acclimol 5l |
| Shelf placement | - |

