



Test case for water analysis in plant aquariums







- Simple and reliable monitoring of the water values in aquariums. Determines: pH value, carbonate hardness, CO2-direct, iron, magnesium freshwater, potassium, phosphate, silicate and nitrate
- Add reagents to water samples, compare with colour charts, read values and look up on protocol sheet if and what action is recommended
- Laboratory comparator system to compensate the inherent water colouring, analysis record sheet, instructions for trouble shooting
- Child-proof reagent bottles, water-proof plastic case, refillable kit
- Package contents: 1 test case for planted aquariums. Incl. 9 reagents, 6 glass vials, 2 syringes, 3 dosing spoons, thermometer, comparator block, 2 plastic cuvettes, colour charts, ball pen, analysis record sheets, instructions



You may also be interested in

You can find a complete overview here: https://www.jbl.de/qr/25511





JBL pH Test 3.0-10.0 Quick test to determine the acidity



JBL pH Test 6.0-7.6 pH test for freshwater aquariums in the range 6.0-7.6



JBL pH Test 7.4-9.0 pH test in aquariums and ponds in the range 7.4-9.0



JBL GH Test Quick test to determine the total hardness











JBL Aquarium Thermometer Float Aquarium thermometer







JBL Testlab ProScape Spare parts







JBL component set for water tests



JBL Testlab empty + inserts

JBL assortment Test **Colour Charts** Colour charts for JBL water



3







Product information

Healthy aquariums with conditions close to nature The right water values are dependent on the fish stock and the plants in the aquarium. Even if the water looks clear it can be contaminated. With bad values diseases or algae can appear in the aquarium. To maintain a healthy aquarium with conditions close to nature it is important to check and adapt the water values regularly. With the JBL Testlab ProScape you can determine the 9 most important water values in your plant aquarium.

Precise measurement of the following values: pH test: acidity of the water from 6.0 to 7.6

CO2 Direct: immediate determination of main plant nutrient carbon dioxide (CO2)

KH test: pH stability of the water (carbonate hardness)

Fe test: iron value to control the fertilisation

NO3 test: monitors the plant nutrient nitrate (another possible reason for algae growth)

PO4 test sensitive: determines the phosphate content (another possible reason for algae growth)

SiO2 test: cause of diatoms

Mg test: plant nutrient magnesium K test: plant nutrient potassium

You can analyse the water values you have measured in the JBL Online Laboratory free of charge and receive recommendations, or have the precise consumption of individual substances in your aquarium individually calculated with the JBL Dosage Calculator. No more following instructions and dosing approximately, in future you can follow the consumption of your aquarium plants and dose precisely!

Further information	
FAQ	~
Blog	~
Press	~
Laboratory/calculator	×
Worth reading	~
Spare parts	~
Video	~
GarantiePlus	×
Instructions	~
QR code	







JBL Testlab ProScape Product details



Article data	
Product name	JBL Testlab ProScape
Art. No.	2551100
EAN number	4014162255112
EAN as barcode	
Content	9 test kits
-	-
Expiry months	24
RRP incl. VAT	-
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	1300 g
Net weight	1300 g
Weight change	1000

Disposal	
Product name	JBL Testlab ProScape
Art. No.	2551100
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	39 g
PPK	-
Plastic small	256 g
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL Testlab ProScape
Art. No.	2551100
Animal species	Arowana, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichlids (South America), Crayfish, Danions, Discus, Dwarf shrimps, Flowerhorn, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Panchaxes, Rainbowfish, Snails, Spiny eels, Tetra, Veiltails, freshwater butterflyfish
Animal size	For all animal sizes



Date: 31.12.2023 Produced by:

Country: Israel

Language: English







Features	
Animal age group	All aquarium fish
Volume of habitat	For all aquariums
Material	-
Food type	-
Colour	multicoloured
Dosage	see instruction manual
Transport conditions	-







Electronic label / illuminant	
Product name	JBL Testlab ProScape
Art. No.	2551100
Ambient temperature	-
Start time	-
Mercury	-
Tube length	-
Service life	-
Lumen	-
CRI value	-
Dimmable	-
Switching cycles	-
PAR value	-
Energy efficiency class	-
UV-A	-
UV-B	-
UV-C	-
Colour temperature	-
Base designation	-

Technical data	
Product name	JBL Testlab ProScape
Art. No.	2551100
Range in litres	-
Range from - to	-
Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	318 mm
Length	382 mm
Width	85 mm
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-
Hose connections suction/in	-
Delivery head	-









Food type	-
Sub product type	-
Dosing	see instruction manual



Date: 31.12.2023 Produced by:









A Safety information

- safety information in accordance with GHS	
Signal word	Danger
Hazard symbol	
Risk phrase	H225: Highly flammable liquid and vapour. H412: Harmful to aquatic life with long lasting effects. H335: May cause respiratory irritation. H341: Suspected of causing genetic defects. H317: May cause an allergic skin reaction. H314: Causes severe skin burns and eye damage. H302: Harmful if swallowed. H261: In contact with water releases flammable gases.
safety note	P101: If medical advice is needed, have product container or label at hand. P102: Keep out of reach of children. P103: Read label before use. P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P301+P310a: IF SWALLOWED: Immediately call a POISON CENTER/ doctor. P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P405: Store locked up. P501a: Dispose of contents/container in accordance with local/regional/ national/international regulations.

Further safety information	
Safety instructions	✓
QR code	



Date: 31.12.2023 Produced by:







Additional information for the specialist trade sector

Article data	
Product name	JBL Testlab ProScape
Art. No.	2551100
VAT	19%
Sales unit (SU)	1
Volume packaging	10.32l
Dimensions (l/w/h)	85 mm/382 mm/318 mm
Layer	24
Pallet	72
Category of products	1
Customs tariff	38220000
Country of origin	DE
Type of packaging	case

PU 1 data	
Product name	JBL Testlab ProScape
Art. No.	2551100
PU 1 material	cardboard gr
PU 1 weight	1700 g
PU 1 lengh	700 mm
PU 1 width/depth	410 mm
PU 1 height	370 mm

PU 2 data	
Product name	JBL Testlab ProScape
Art. No.	2551100
PU 2 material	None
PU 2 weight	-
PU 2 lengh	-
PU 2 width/depth	-
PU 2 height	-

Trade data	
Product name	JBL Testlab ProScape
Art. No.	2551100
Till receipt text	Testlab ProScap
Shelf placement	-



31.12.2023 Date: Produced by:

10