



Main food for tortoises 10 - 50 cm in size

Suitable for:









- Especially for the nutritional requirements of tortoises: food sticks made up of a mixture of purely plantbased raw materials
- To maintain the animals' health: selected herbs and plants
- Ideal for the digestive system: high fibre content, low protein content
- Healthy growth without shell problems, stabilises and strengthens immunity thanks to vitamin complex and vitamin C
- Contents: staple food for tortoises, Agivert. Can be stored for 3 years unopened, use within 3 months after opening



You may also be interested in

You can find a complete overview here: https://www.jbl.de/qr/70331





JBL Tortoise Sun Terra Vitamins for tortoises



JBL TerraVit fluid Vitamins and trace elements for terrarium animals



JBL MicroCalcium Mineral supplementary food for all reptiles



JBL TerraVit Powder Vitamins and trace elements for terrarium animals







JBL Agivert Accessories





JBL Tortoise Sun Terra Vitamins for tortoises



JBL TerraVit fluid Vitamins and trace elements for terrarium animals



JBL MicroCalcium Mineral supplementary food for all reptiles



JBL TerraVit Powder Vitamins and trace elements for terrarium animals




JBL Tortoise shell care Shell care for tortoises











Especially for the nutritional requirements of tortoises
To remain healthy animals which are predominantly herbivore need food which is low in protein and rich in fibre. The food for tortoises deliberately contains only vegetable ingredients, high in fibre and low in proteins.

Feeding recommendations Several times a day feed as much as can be eaten in a short time. For adult animals 4 to 5 feeds per week are sufficient. Feed directly, soaked in small quantities of water or with fresh meadow herbs.

| Further information | |
|-----------------------|----------|
| FAQ | ~ |
| Blog | ~ |
| Press | ✓ |
| Laboratory/calculator | × |
| Worth reading | ~ |
| Spare parts | × |
| Video | × |
| GarantiePlus | × |
| Instructions | ~ |
| QR code | |




31.12.2023 Date: Produced by:







JBL Agivert Product details



| Article data | | | |
|--------------------------|--------------------|--------------------|-----------------|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 I |
| Art. No. | 7033100 | 7033200 | 7033300 |
| EAN number | 4014162703316 | 4014162703323 | 4014162703330 |
| EAN as barcode | | | |
| Weight | 42 g | 105 g | 420 g |
| Size | 100 ml | 250 ml | 1000 ml |
| Expiry months | 36 | 36 | 36 |
| RRP incl. VAT | 4,09 € | 6,31 € | 18,44 € |
| Base price | 95.12 € | 60.1 € | 43.9 € |
| Nominal filling quantity | 43 g | 105 g | 420 g |
| Base quantity | 1 kg | 1 kg | 1 kg |
| Gross weight | 58.6 g | 132.3 g | 498.9 g |
| Net weight | 42 g | 105 g | 420 g |
| Weight change | 1000 | 1000 | 1000 |

| Disposal | | | | |
|-------------------------|--------------------|--------------------|-----------------|--|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 l | |
| Art. No. | 7033100 | 7033200 | 7033300 | |
| Green dot | ✓ | ✓ | ✓ | |
| Group electronic waste | - | - | - | |
| Disposal weight | - | - | - | |
| Battery type | - | - | - | |
| Battery return | - | - | - | |
| Battery rechargeable | - | - | - | |
| Disposal weight battery | - | - | - | |
| Non-returnable glass | - | - | - | |
| PPK | - | - | - | |
| Plastic small | 16.6 g | 27.3 g | 78.9 g | |
| Plastic large | - | - | - | |
| Disposal weight metal | 0 g | 0 g | 0 g | |

| Features | | | | |
|-------------------|---|---|---|--|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 l | |
| Art. No. | 7033100 | 7033200 | 7033300 | |
| Animal species | European tortoises, Leopard tortoises, Red-footed tortoises, Tropical tortoises | European tortoises, Leopard tortoises, Red-footed tortoises, Tropical tortoises | European tortoises, Leopard tortoises, Red-footed tortoises, Tropical tortoises | |
| Animal size | 10 - 50 cm | 10 - 50 cm | 10 - 50 cm | |
| Animal age group | Medium-sized tortoises | Medium-sized tortoises | Medium-sized tortoises | |
| Volume of habitat | For all terrariums | For all terrariums | For all terrariums | |



31.12.2023 Date: Produced by:







| Features | | | | |
|-----------------------|---|--|--|--|
| Material | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition | You will find details of the materials used under Details & Composition | |
| Food type food sticks | | food sticks | food sticks | |
| Colour | green | green | green | |
| Dosage | Feed young, growing turtles in the water area, as much JBL Calcil as can be eaten in about 10 minutes, 1 - 2 times a day, depending on size., 2 - 3 feeds per week in the same way are sufficient for adult turtles | Feed young, growing turtles in the water area, as much JBL Calcil as can be eaten in about 10 minutes, 1 - 2 times a day, depending on size. 2 - 3 feeds per week in the same way are sufficient for adult turtles | Feed young, growing turtles in the water area, as much JBL Calcil as can be eaten in about 10 minutes, 1 - 2 times a day, depending on size. , 2 - 3 feeds per week in the same way are sufficient for adult turtles | |
| Transport conditions | - | - | - | |







| Electronic label / illuminant | | | | |
|-------------------------------|--------------------|--------------------|-----------------|--|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 I | |
| Art. No. | 7033100 | 7033200 | 7033300 | |
| Ambient temperature | - | - | - | |
| Start time | - | - | - | |
| Mercury | - | - | - | |
| Tube length | - | - | - | |
| Service life | - | - | - | |
| Lumen | - | - | - | |
| CRI value | - | - | - | |
| Dimmable | - | - | - | |
| Switching cycles | - | - | - | |
| PAR value | - | - | - | |
| Energy efficiency class | - | - | - | |
| UV-A | - | - | - | |
| UV-B | - | - | - | |
| UV-C | - | - | - | |
| Colour temperature | - | - | - | |
| Base designation | - | - | - | |

| Technical data | | | |
|-------------------------------|--------------------|--------------------|-----------------|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 l |
| Art. No. | 7033100 | 7033200 | 7033300 |
| Range in litres | - | - | - |
| Range from - to | - | - | - |
| Range in days | - | - | - |
| Range tank length | - | - | - |
| Output in watts | - | - | - |
| Output per hour | - | - | - |
| Output per day | - | - | - |
| Height | - | 120 mm | - |
| Length | - | 60 mm | - |
| Width | - | 60 mm | - |
| Diameter | - | 60 mm | - |
| Voltage | - | - | - |
| For | - | - | - |
| T8 26mm (watt) | - | - | - |
| T5 16mm (watt) | - | - | - |
| Size | - | - | - |
| Content for | - | - | - |
| Filter container volume | - | - | - |
| Volume filter media | - | - | - |
| Hose connections pressure/out | - | - | - |
| Hose connections suction/in | - | - | - |
| Delivery head | - | - | - |



Country: Moldova

Language: English





| Food type | food sticks |
|------------------|---|
| Sub product type | complete food for tortoises |
| Dosing | Feed young, growing turtles in the water area, as much JBL Calcil as can be eaten in about 10 minutes, 1 - 2 times a day, depending on size. ,2 - 3 feeds per week in the same way are sufficient for adult turtles |


Analytical components



| | Crude protein | 12.5 % | |
|--|---------------|--------|--|
| | Fat content | 2.5 % | |
| | Crude fibre | 22 % | |
| | Crude ash | 8.5 % | |

III Composition





I ★ Additives

| Antioxidants | |
|------------------------------------|--|
| E 306 (natural vitamin E extracts) | |

| Vitamins, provitamins and other chemically defined substances having a similar effect (per 1000 g) | | | |
|--|--|--|--|
| Vitamin A 8000 I. E. | | | |
| Vitamin D₃ 800 I. E. | | | |
| Vitamin E 140 mg | | | |
| Vitamin C(stable) 100 mg | | | |











Additional information for the specialist trade sector

| Article data | | | |
|----------------------|--------------------|--------------------|--------------------|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 I |
| Art. No. | 7033100 | 7033200 | 7033300 |
| VAT | 7% | 7% | 7% |
| Sales unit (SU) | 6 | 6 | 3 |
| Volume packaging | 0.211 | 0.43l | 1.72l |
| Dimensions (l/w/h) | 50 mm/50 mm/85 mm | 60 mm/60 mm/120 mm | 98 mm/98 mm/179 mm |
| Layer | 672 | 648 | 192 |
| Pallet | 2688 | 1944 | 576 |
| Category of products | 1 | 1 | 1 |
| Customs tariff | 12149090 | 12149090 | 12149090 |
| Country of origin | DE | DE | DE |
| Type of packaging | can | can | can |

| PU 1 data | | | | |
|------------------|--------------------|--------------------|-----------------|--|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 l | |
| Art. No. | 7033100 | 7033200 | 7033300 | |
| PU 1 material | film gr | film gr | film gr | |
| PU 1 weight | 2.4 g | 4.5 g | 9 g | |
| PU 1 lengh | 92 mm | 117 mm | 103 mm | |
| PU 1 width/depth | 142 mm | 174 mm | 300 mm | |
| PU 1 height | 85 mm | 120 mm | 179 mm | |

| PU 2 data | | | | |
|------------------|--------------------|--------------------|-----------------|--|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 l | |
| Art. No. | 7033100 | 7033200 | 7033300 | |
| PU 2 material | cardboard gr | cardboard gr | cardboard gr | |
| PU 2 weight | 357 g | 430 g | 510 g | |
| PU 2 lengh | 300 mm | 255 mm | 315 mm | |
| PU 2 width/depth | 395 mm | 390 mm | 415 mm | |
| PU 2 height | 180 mm | 370 mm | 340 mm | |

| Trade data | | | | | |
|-------------------|--------------------|--------------------|-----------------|--|--|
| Product name | JBL Agivert 100 ml | JBL Agivert 250 ml | JBL Agivert 1 l | | |
| Art. No. | 7033100 | 7033200 | 7033300 | | |
| Till receipt text | Agivert 100ml | Agivert 250ml | Agivert 1I | | |
| Shelf placement | - | - | - | | |



31.12.2023 Date: Produced by: