

Spare filter foam for Cristal Profi m module

Suitable for: < 🖈 🐋





- Clean, healthy water: spare filter foam for the extension module of the des JBL CristalProfi m internal filter
- Exact fit
- For module 13 x 9 cm.
- Contains two spare filter foams
- Package contents: spare filter foams for JBL Cristal Profi m module, CristalProfi m greenline modul filterpad



You may also be interested in You can find a complete overview here: https://www.jbl.de/qr/60968





JBL CristalProfi m greenline Module Filter module to extend CristalProfi m





Product information

Quality of foam:

Foams used in aquaristics with live animals must always be subjected to a pollutant test, as most plastics release problematic plasticisers into the water. All the plastics used by JBL have passed this test and do NOT release any harmful substances into the water.

Cleaning of plastics:

For many months it is sufficient to rinse the dirty filter sponges under lukewarm running water and then inoculate them with filter bacteria (JBL FilterStart) before using them again. However, as mechanical cleaning NEVER leads to complete cleaning and dirt will remain in the filter sponge, it is recommended you replace the filter sponges about every 6 months.

you replace the filter sponges about every 6 months. If you place a freshly cleaned filter sponge in a glass of clean water overnight, you will be able to smell how much dirt is still present in the supposedly clean sponge the next morning!

Inoculating filter material with cleansing bacteria:

Initial colonisation of filter material with cleansing bacteria: Initial colonisation of filter material with cleansing bacteria takes about 1-2 weeks in freshwater and twice as long in seawater. If the filter material was already in use, the time for a "recolonisation" is doubled, as recolonisation does not take place so easily on dead bacterial lawns. The addition of concentrated cleaning bacteria such as in JBL FilterStart shortens the biological start-up time to about 24 hours.

Another reason is to fill the existing settlement space in the filter material with beneficial bacteria so that pathogenic germs find no place to spread. This has now been scientifically proven.

Further information	
FAQ	~
Blog	~
Press	~
Laboratory/calculator	×
Worth reading	~
Spare parts	×
Video	~
GarantiePlus	×
Instructions	~
QR code	



Article data	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800
EAN number	4014162609687
EAN as barcode	
Content	2 piece(s)
-	-
Expiry months	-
RRP incl. VAT	7,29€
Base price	-
Nominal filling quantity	-
Base quantity	1
Gross weight	56.5 g
Net weight	14 g
Weight change	1000

Disposal	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800
Green dot	✓
Group electronic waste	-
Disposal weight	-
Battery type	-
Battery return	-
Battery rechargeable	-
Disposal weight battery	-
Non-returnable glass	-
РРК	28 g
Plastic small	-
Plastic large	-
Disposal weight metal	0 g

Features	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800

Country: Thailand

Language: English

TECHNICAL PRODUCTS TECHNIK TECHNIQUE



Features	
Animal species	Armored catfish, Arowana, Artemia, Axolotl, Barbels, Bettas, Bichirs/reedfish, Blowfish, Catfish, Cichilds (South America), Clawed frogs, Clawed shrimps, Corals, Crayfish, Crustaceans, Danions, Discus, Dwarf crayfish, Dwarf shrimps, Flowerhorn, Gill maggots, Gobies, Goldfish, Gouramis, Guppy, Juvenile fish, Killifish, Livebearers, Loaches, Mussels, Newts, Panchaxes, Rainbowfish, Snails, Spiny eels, Terrapins, Tetra, Tropical terrapins, Veiltails, blood parrot cichlids, freshwater butterflyfish, turtle
Animal size	1 - 5 cm
Animal age group	Small aquarium fish
Volume of habitat	20 - 80 (20 - 100 cm)
Material	plastic (ABS) / polyether
Food type	-
Colour	black
Dosage	-
Transport conditions	-





JBL CP m greenline FilterPad module

Electronic label / illuminant		
Product name	JBL CP m greenline FilterPad module	
Art. No.	6096800	
Ambient temperature	-	
Start time	-	
Mercury	-	
Tube length	-	
Service life	-	
Lumen	-	
CRI value	-	
Dimmable	-	
Switching cycles	-	
PAR value	-	
Energy efficiency class	-	
UV-A	-	
UV-B	-	
UV-C	-	
Colour temperature	-	
Base designation	-	
Technical data		
Product name	JBL CP m greenline FilterPad module	
Art. No.	6096800	
Range in litres	-	
Range from - to	-	
Range in days	-	
Range tank length	-	
Output in watts	-	

Range in days	-
Range tank length	-
Output in watts	-
Output per hour	-
Output per day	-
Height	-
Length	-
Width	-
Diameter	-
Voltage	-
For	-
T8 26mm (watt)	-
T5 16mm (watt)	-
Size	-
Content for	-
Filter container volume	-
Volume filter media	-
Hose connections pressure/out	-
Hose connections suction/in	-



TECHNICAL PRODUCTS TECHNIK TECHNIQUE





Technical data	
Delivery head	-





JBL CP m greenline FilterPad module

Food type	-
Sub product type	-
Dosing	-





Additional information for the specialist trade sector

Article data	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800
VAT	19%
Sales unit (SU)	3
Volume packaging	0.921
Dimensions (l/w/h)	45 mm/150 mm/137 mm
Layer	378
Pallet	756
Category of products	3
Customs tariff	84219990
Country of origin	CN
Type of packaging	Faltschachtel/Karton

PU 1 data	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800
PU 1 material	film gr
PU 1 weight	1 g
PU 1 lengh	100 mm
PU 1 width/depth	135 mm
PU 1 height	140 mm

PU 2 data	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800
PU 2 material	cardboard gr
PU 2 weight	1050 g
PU 2 lengh	730 mm
PU 2 width/depth	320 mm
PU 2 height	440 mm

Trade data	
Product name	JBL CP m greenline FilterPad module
Art. No.	6096800
Till receipt text	Profi m Mod Pad
Shelf placement	-

